

# Sono-Eyes: Get Sonographic Vision

Judith Gull, Samuel Hertig, Christoph Krautz, Andreas Pedroni,  
Ueli Peter, Alessandro Pianezzi, Till Bay

Incremed AG, Zürich

## Advantages and Challenges of Ultrasound Imaging

Ultrasound offers precise and safe real-time imaging at low cost. Point-of-care ultrasound enables clinicians to treat patients wherever they are located (Ref. 1). To fully harness the capabilities of ultrasound, the following challenges must be solved:

- Overcome the disconnect between the ultrasound image on screen and the location of the probe
- Get rid of screens that are bulky and need sterilization in operating rooms
- Provide real time needle guidance for ultrasound procedures such as needle insertions (Ref. 2)



Figure 1. Microsoft HoloLens mixed reality headset.

## Sono-Eyes Facilitates Needle Insertions

Sono-Eyes provides measurable value for ultrasound training (Ref. 3):

- A preliminary study demonstrates faster needle insertions for non-radiologists (Fig. 5A) when using Sono-Eyes
- Non-radiologists also have a higher percentage of single-pass insertions with Sono-Eyes compared to using classical US (Fig. 5B)

## Outlook

- Integrate ultrasound devices fully such that all device functions can be operated through the mixed reality of Sono-Eyes
- Planned *in-vivo* studies
- Advancements in headsets will further accelerate the integration of mixed reality for daily clinical use




Figure 2. Sono-Eyes *in-situ* ultrasound.

## Sono-Eyes: *In-situ* Mixed Reality Ultrasound

- Sono-Eyes is a ultrasound solution designed to eliminate current shortcomings with the help of mixed reality. The software runs on a headset worn by the clinician (Fig. 1)
- The clinician can now see the ultrasound image at its anatomical position, which completely eliminates the disconnect between screen, patient and clinician (Fig. 2)
- A virtual screen can be enabled at any time, which provides a larger image that follows the clinician's gaze without obstructing the line of sight to the patient
- Operating an ultrasound device hands-free is key when performing sterile interventions: Settings are adjusted using gazable buttons (Fig. 3)
- To assist ultrasound-guided needle insertions, Sono-Eyes provides visual guidance for an optimal needle path (Fig. 4)



Figure 3. Gazable button.




Figure 4. A) Frontal view of needle guides. B) Lateral view.




Figure 5. A) With non-specialists, Sono-Eyes reduces needle insertion time for ultrasound-guided procedures when compared to using a traditional ultrasound screen. B) Sono-Eyes increases the percentage of single-pass insertion attempts.

## References and Acknowledgments

- 1) Moore, CL & Copel, JA; "Point-of-Care Ultrasonography", N Engl J Med, 2011
- 2) Scholten, HJ et al; "Improving needle tip identification during ultrasound-guided procedures in anaesthetic practice", Anaesthesia, 2017
- 3) Farshad-Amacker NA, Bay T, Roskopf AB, Spirig JM, Wanivenhaus F, Pfirrmann CWA, Farshad M, manuscript in preparation

The authors would like to thank I. Kanpolat and R. Schiek for assistance with photography, as well as N. Farshad-Amacker and M. Farshad for helpful discussions. Correspondence to [till.bay@incremed.com](mailto:till.bay@incremed.com)

Der **Balgrist**

[www.incremed.com](http://www.incremed.com)

